

**Author:** Department of Medical Sciences, Ministry of Public Health, Thailand

**Technique:** Real Time RT-PCR

**Methodology:** The protocol is designed in order to target the *N* gene to detection of 2019-nCoV infection.

**Material:**

- NucleoSpin RNA Virus, Mini kit for viral RNA from cell-free fluids (Macherey-Nagel; Cat. No 740956)
- Superscript™ III Platinum One-Step qRT-PCR Kit (Invitrogen™; Cat No. 11732020 or 11732088)
- Primers and Probe sequences

WH-NIC N-F: 5'-CGTTTGGTGGACCCTCAGAT-3'

WH-NIC N-R: 5'-CCCCACTGCGTTCTCCATT-3'

WH-NIC N-P: 5'-FAM-CAACTGGCAGTAACCA-BQH1-3'

**Protocol:**

1. Viral RNA extraction, performed according to NucleoSpin RNA Virus manufacturer's instructions
2. Preparation of one-step multiplex RT-PCR master mix as below. The reactions should be set up on ice. Negative and positive controls should be performed.  
[Note: Primers and probe are used at different concentrations]

Reagent	Volume (µL) per reaction
H <sub>2</sub> O (RNase free)	5.5
2x PCR Reaction Mix	12.5
SuperScript-Taq Mix	0.5
Primer WH-NIC N-F (40µM)	0.5
Primer WH-NIC N-R (40µM)	0.5
Probe WH-NIC N-P (10µM)	0.5
Total reaction mix	20
Template RNA, add	5
Total volume	25

3. Set the follow RT-PCR conditions:

Temperature (°C)	Time	No of Cycles
50	30 min	1
95	2 min	
95	15 sec	45
55	30 sec	

**Results Evaluation:** Negative control should be undetected. A positive control should be detected with  $Ct \leq 38$ .