

**Author:** Institute Pasteur, Paris

**Technique:** Real Time RT-PCR

**Methodology:** The protocol is designed in order to detect SARS-CoV-2 with two *RdRp* gene targets (IP2 and IP4) in a multiplex reaction. A singleplex confirmatory assay analysis of *E* gene is also performed.

**Material:**

- Kit Extraction NucleoSpin Dx Virus (Macherey Nagel; 740895.50)
- SuperScript™ III Platinum® One-Step Quantitative RT-PCR System (Invitrogen; 1732-020)
- LightCycler 480 (96)
- Positive control – *in vitro* transcribed RNA derived from strain BetaCoV\_Wuhan\_WIV04\_2019 (EPI\_ILS\_402124) (available upon request; grippe@pasteur.fr)

- Primers and Probes sequences

*RdRp* gene / *nCoV\_IP2* (108 bp)

nCoV\_IP2-12669Fw: 5'-ATGAGCTTAGTCCTGTTG-3'

nCoV\_IP2-12759Rv: 5'-CTCCCTTTGTTGTGTTGT-3'

nCoV\_IP2-12696bProbe(+): 5'-HEX-AGATGTCTGTGCTGCCGGTA-BHQ-1-3'

*RdRp* gene / *nCoV\_IP4* (107 bp)

nCoV\_IP4-14059Fw: 5'-GGTAACTGGTATGATTTTCG-3'

nCoV\_IP4-14146Rv: 5'-CTGGTCAAGGTTAATATAGG-3'

nCoV\_IP4-14084Probe(+): 5'-FAM-TCATACAAACCACGCCAGG-BHQ-1-3'

*E* gene / *E\_Sarbeco* (125 bp)

E\_Sarbeco\_F1: 5'-ACAGGTACGTTAATAGTTAATAGCGT-3'

E\_Sarbeco\_R2: 5'-ATATTGCAGCAGTACGCACACA-3'

E\_Sarbeco\_P1: 5'-FAM-ACACTAGCCATCCTTACTGCGCTTCG-BBQ-1-3'

**Protocol:**

1. Viral RNA extraction, performed according to NucleoSpin Dx Virus kit manufacturer's instructions. RNA extracted from 100 µL of original sample, is eluted in 100 µL of elution buffer.

2. Preparation of real-time RT-PCR master mix as presented below. While the analysis of *E* gene should be performed in a simplex mix, the primer sets of *nCoV\_IP2* and *nCoV\_IP4* can be multiplexed. Negative and positive controls should be performed.

[Other kits can be used to perform RT-PCR. For the matter, the optimized concentrations of the reagents are described below for each reaction mixture]

Simplex Mix		
Reagent	Volume ( $\mu$ L) per reaction	Final Concentration
H <sub>2</sub> O PPI	3.60	
Reaction Mix 2x	12.50	3 mM Mg
MgSO <sub>4</sub> (50mM)	0.40	0.8 mM Mg
SuperScript III/Platinum Taq Mix	1.00	
Primer FW (10 $\mu$ M)	1.00	0.4 $\mu$ M
Primer RV (10 $\mu$ M)	1.00	0.4 $\mu$ M
Probe (10 $\mu$ M)	0.50	0.2 $\mu$ M
Total reaction mix	20	
Template RNA, add	5	
Total volume	25	

Multiplex Mix		
Reagent	Volume ( $\mu$ L) per reaction	Final Concentration
H <sub>2</sub> O PPI	1.30	
Reaction Mix 2x	12.50	3 mM Mg
MgSO <sub>4</sub> (50mM)	0.40	0.8 mM Mg
SuperScript III/Platinum Taq Mix	1.00	
Primer FW (10 $\mu$ M)	1.00	0.4 $\mu$ M
Primer RV (10 $\mu$ M)	1.00	0.4 $\mu$ M
Primer FW (10 $\mu$ M)	1.00	0.4 $\mu$ M
Primer RV (10 $\mu$ M)	1.00	0.4 $\mu$ M
Probe (10 $\mu$ M)	0.4	0.16 $\mu$ M
Probe (10 $\mu$ M)	0.4	0.16 $\mu$ M
Total reaction mix	20	
Template RNA, add	5	
Total volume	25	

3. Set the follow RT-PCR conditions:

Temperature (°C)	Time	No of Cycles
55	20 min	1
95	3 min	
95	15 sec	50
58	30 sec	
40	30 sec	1

**Results Evaluation:** Ct values for the reactions are presented below. The values may vary from instrument to instrument by up to 2 cycles, while the interval between two dilutions steps is constant ( $\Delta$ Ct).

Ct values for the performed reactions			
	Multiplex		Singleplex
RNA copies of transcript	nCoV_IP2	nCoV_IP4	E_Sarbeco
1,00E+07	21.67	21.97	24.72
1,00E+06	24.97	25.12	28.19
1,00E+05	28.00	27.88	30.96
1,00E+04	31.84	30.51	33.33